

PCR protocol

RBRC 02436

Strain name B6.Cg-Tg(CAG-Cadm1)

Primer	primer_name	primer_seq	length		
1	Sgi-center-R	gtccgaccactcctccacctctgattggcc	30		
2	IGSF4-86F	tccggctcctgctgttgccttt	24		
3	oIMR0042	CTAGGCCACAGAATTGAAAGATCT	24	Internal control (IL-2)	JAX Genotyping Protocol
4	oIMR0043	GTAGGTGGAAATTCTAGCATCATCC	25	Internal control (IL-2)	JAX Genotyping Protocol

Reaction Components

Enzyme TaKaRa LA Taq with GC buffer (Cat. No. RR02AG)

Reaction Conditions

	Conc	Run 1	Run 2	Run 3	Run 4
H2O		2.8			
PCR buffer	2 X GC I	10.0			
dNTP	2.5 mM	2.0			
primer 1	10 pmol/ μ l	1.0			
primer 2	10 pmol/ μ l	1.0			
primer 3	10 pmol/ μ l	1.0			
primer 4	10 pmol/ μ l	1.0			
DNA		1.0			
Taq		0.2			
total		20			

step	temp (°C)	time (sec)	note
1	94	180	
2	94	30*	
3	55	30*	
4	72	60*	*Repeat for 30 cycles
5	72	300	

Marker

100 bp ladder (Invitrogen, Cat. 15628-050)

Gel

2% agarose in 1 X TAE (100 V, 30 min)

PCR products

primer_set	product_size	product
primer 1	521 bp	Tg
primer 3	324 bp	Internal control (IL-2)

Comment